C612M Digital Torque Tester

The application torque and removal torque of caps or closures of bottles, spout bags and flexible tube packages are important online and offline production parameters for the manufacturers, they may influence the transportation and consumption of the products. C612M Digital Torque Tester is the first generation of automatic torque tester with manipulator, which is one professional testing instrument that can be used to measure the application torque and removal torque of the caps of bottles, spout bags and flexible tube packages. C612M, featuring great stability and high precision, is an indispensable testing instrument in manufacturing process.

Product Features^{Note1}

Dual Test Modes & Automatic Test with Manipulator

- Dual test modes including application torque test and removal torque test
- Innovative patented manipulator technique supports automatic gripping, opening and fastening, which minimizes the error caused by manual operations and improves the testing accuracy and repeatability
- Gripping force, fastening force and turning speed of the cap can be adjusted freely
- Automatic gripping manipulator with deviation less than 0.01Nm, much better than manual fastening
- Intelligent designs of overload protection, automatic zeroing and error alert for safe test operation
- Manual test or automatic test can be selected

High Testing Accuracy & Low Testing Lower Limit

- Specimens with torque value lower than 0.005Nm is testable. The testing resolution is as high as 0.0001Nm, providing high accuracy and good repeatability.
- The peak value can be maintained automatically, which ensures the recording of the accurate test results
- Automatic judgement of peak value meets the testing requirements for various specimens
- Components and parts are supplied by world well-known manufacturers, which guarantee the stable performance
- Imported pneumatic control system with ultra-long service life and ultra-low failure rate ensures the testing accuracy

New Patented Intelligent Touch Operating System

- Industrial grade touch screen, one-button operation, user-friendly interface, remote upgrade and maintenance
- Chinese and English languages can be selected
- Test curves are displayed in real time. Intelligent data statistics supports quick review of test results.
- Automatic data storage and power failure memory functions
- Historical data can be reviewed and printed easily
- Storage up to 1200 test records
- Test data is expressed in eight international units for reference and comparison
- Multi-level user privilege management, log in with password

 Labthink Instruments Co., Ltd.
 144 Wuyingshan Road, Jinan, P.R.China (250031)
 Phone: +86-531-85068566
 FAX: +86-531-85062108

 Labthink International, Inc.
 200 River's Edge Drive, Medford, MA, 02155, U.S.A.
 Phone: +1-617-830-2190
 FAX: +1-781-219-3638



Labthink®



- Micro-printer and USB data ports for data exporting and data transfer (optional)
- The instrument conforms to the requirements of GMP (optional)
- Labthink's unique DataShiledTM System for data management and connecting with information system (optional)

Test Standards^{Note1}

GB/T 17876, ASTM D2063, ASTM D3198, ASTM D3474, BB/T 0025, BB/T 0034

Applications^{Note1}

Derie	Bottle Packages	Application and removal torque test of the bottle caps for bottled packages for food, pharmaceuticals, e.g. beverage or tablet bottles
Basic Applications	Flexible Tube Packages	Application and removal torque test of flexible tube packages for food, pharmaceuticals and cosmetics, e.g. eye-drop, hand cream and shoe polish
Extended Applications	Screws	Application and removal torque test of the screws
	Vacuum Flasks and	Application and removal torque test of the caps of vacuum flasks and
	Vacuum Cups	vacuum cups

Technical Specifications^{Note2}

Specifications	C612M	
	5 Nm (Standard)	
Load Cell Capacity	20 Nm (Optional)	
	40 Nm (Optional)	
A	Indicated Value ±0.5% (10%-100% of Load cell capacity)	
Accuracy	$\pm 0.05\%$ FS (0%-10% of Load cell capacity)	
Resolution	0.0001 Nm	
Body Clamp Range	Φ5 mm ~ Φ170 mm	
Cap Clamp Range	$\Phi 10 \text{ mm} \sim \Phi 80 \text{ mm}$	
Bottle Height	20 mm ~ 400 mm	
Specimen Clamping and	Pneumatic	
Turning Way		
Maximum Auto Removal	2 New (constantion considering of the set of	
/ Application Torque	2 Nm (customization available for other specs)	
Gas Supply	Air (Not in supply scope)	
Gas Supply Pressure	0.7 Mpa (101.5psi)	
Statistical Analysis	0000	
Allowance	0 ~ 999	
Instrument Dimension	550 mm (L) x 365 mm (W) x 1150 mm (H)	

 Labthink Instruments Co., Ltd.
 144 Wuyingshan Road, Jinan, P.R.China (250031)
 Phone: +86-531-85068566
 FAX: +86-531-85062108

 Labthink International, Inc.
 200 River's Edge Drive, Medford, MA, 02155, U.S.A.
 Phone: +1-617-830-2190
 FAX: +1-781-219-3638

Labthink®

Power Supply	220VAC±10% 50Hz / 120VAC±10% 60Hz
Net Weight	39 kg

Configurations

Standard Configurations	Instrument, Gripping Post (4pcs), Gripping Block (1 pair), Calibration Kit (not
Standard Configurations	include weight), Φ 4mm PU Tubing (2m)
Ontional Porta	Micro-printer, Professional Software, Air Compressor, GMP System,
Optional Parts	DataShield ^{TM Note3}
Nata	1. The gas supply port of the instrument is $\Phi 4 \text{ mm PU Tubing}$;
Note	2. Customers need to prepare gas supply.

Note 1: The described test standards, applications and product features should be in line with Technical Specifications.

Note 2: The parameters in the table are measured by professional operators in Labthink laboratory under strictly controlled laboratory conditions.

Note 3: DataShieldTM provides safe and reliable data application support. Multiple Labthink instruments can share one single DataShieldTM system which can be configured as required.

Please Note: Labthink is always dedicated to the innovation and improvement of product performance and function. Therefore, technical specifications are subject to change without further notice. Please visit our website at www.labthink.com for the latest updates. Labthink reserves the rights of final interpretation and revision.

Labthink Instruments Co., Ltd. 144 Wuyingshan Road, Jinan, P.R.China (250031) Phone: +86-531-85068566 FAX: +86-531-85062108 Labthink International, Inc. 200 River's Edge Drive, Medford, MA, 02155, U.S.A. Phone: +1-617-830-2190 FAX: +1-781-219-3638